


C2075-0517

Determination of Tricyclic Antidepressants via LC-UV using ChraSil B100 C18, 150mm x 4.6mm, 5µm

Column	: ChraSil B100 C18, 150mm x 4.6mm, 5µm	Peak Identification 1. Nordoxepin 2. Doxepin 3. Desipramine 4. Imipramine 5. Nortriptyline 6. Amitriptyline 7. Trimipramine 8. Clomipramine
Part No.	: C1505-B5032	
Mobile Phase	: A: H ₂ O + 0.1% HCOOH, B: CH ₃ CN + 0.1% HCOOH	
Gradient	: 30% - 40% B in 15 minutes	
Flow Rate	: 1.0 mL/min	
Temperature	: Ambient	
Detection	: UV, 254 nm	
 image 1627371932648.png		
ChraNotes at https://chranote.com/C2075-0517		

Revision #8

Created 27 July 2021 07:36:29 by Admin

Updated 6 October 2021 06:28:13 by Admin