

C2086-0508

Determination of Antihistamines by HPLC-UV together with ChraSil D120 C8, 150mm x 4.6mm, 5µm

Column	: ChraSil D120 C8, 150mm x 4.6mm, 5µm	Peak Identifi
Part No.	: C1440-D5032	1. Pheniram
Mobile Phase	: CH ₃ OH:5 mM NH ₄ HCO ₃ , pH10 - 75:25	2. Doxylami
Flow Rate	: 1.0 mL/min	3. Pyrilamin
Temperature	: Ambient	4. Phenothia
Detection	: UV, 254 nm	5. Diphenhy
		6. Loratadin
		7. Prometha
		8. Chlocylizi

image-1627384964558.png

ChraNotes

Revision #3

Created 27 July 2021 11:18:17 by Admin

Updated 27 July 2021 11:29:33 by Admin