

C2091-0513

Determination of Phenoxyacid Herbicides via LC-UV by ChraSil D120 C8,150mm x 4.6mm, 5µm

Column : ChraSil D120 C8,150mm x 4.6mm, 5µm

Part No. : [C1440-D5032](#)

Mobile Phase : CH3CN: H2O + 0.1% HCOOH - 50:50

Flow Rate : 1.0 mL/min

Temperature : Ambient

Detection : UV, 220 nm

Peak Identification

1. Phenoxyac
2. o-Chloroph
3. p-Chloroph
4. 2,3-Dichlor
5. 2,4-Dichlor
6. 2,4,5-Trich
7. 2,4,5-Trich

[image-1627386106132.png](#)

ChraNotes

Revision #1

Created 27 July 2021 11:40:26 by Admin

Updated 27 July 2021 11:42:38 by Admin