

# C1660-D5034

|                    |  |
|--------------------|--|
| Description        | ChraSil HPLC Column  |
| Manufacturer       | Chratec  |
| Internal Diameter  | 4.6  |
| Column Length      | 250  |
| Particle Size      | 5.0  |
| Modification       | ODS-2  |
| Pressure Stability | 400 bar  |
| Base Material      | Silica, spherical, porous  |
| Packing brand      | ChraSil  |
| USP L Code         | L1   |
| Part Number        | C1660-D5034  |
| pH Range#          | 2 - 8  |
| Carbon Load#       | - (See <a href="https://chratesolutions.com">https://chratesolutions.com</a> ) |
| Endcap             | #  |
| Pore Size          | -  |
| Surface Area(m2/g) | -  |
| Column Group#      | -  |
| ChraSil URL#       | <a href="#">C1660-D5034</a> at <a href="#">https://chratesolutions.com</a>     |
| Installation       | Please see instrument manual   |
| Replacement        | Beyond routine maintenance   |

|              |   |
|--------------|---|
| Care and Use | <a href="https://chrates.com">https://chrates.com</a> |
| Warranty#    | See ChraContact at                                    |

---

Revision #1

Created 30 July 2021 05:49:26 by Admin

Updated 30 July 2021 05:49:33 by Admin